

A complex network diagram with numerous nodes and connecting lines, resembling a molecular structure or a data network, serves as the background for the entire page.

# Join Application Operation Manual

Android


# Contents

---

1.	<a href="#">Installation &amp; Registration Process</a>	• • • • • 3
2.	<a href="#">Caution for Using</a>	• • • • • 6
3.	<a href="#">Change registered information</a>	• • • • • 7
4.	<a href="#">Join Reinstallation</a>	• • • • • 9
5.	<a href="#">When you forgot password</a>	• • • • • 10
6.	<a href="#">Contacts &amp; Groups</a>	• • • • • 11
7.	<a href="#">Direct Chat</a>	• • • • • 13
8.	<a href="#">Direct Chat Function ①</a>	• • • • • 14
9.	<a href="#">Direct Chat Function ②</a>	• • • • • 15
10.	<a href="#">Direct Chat Function ③</a>	• • • • • 16
11.	<a href="#">Group Chat</a>	• • • • • 17
12.	<a href="#">Group Chat Function ①</a>	• • • • • 18
13.	<a href="#">Group Chat Function ②</a>	• • • • • 19
14.	<a href="#">Group Chat Function ③</a>	• • • • • 20
15.	<a href="#">Notifications ①</a>	• • • • • 21
16.	<a href="#">Notifications ②</a>	• • • • • 22
17.	<a href="#">DICOM Image Acquisition</a>	• • • • • 23
18.	<a href="#">Belonging to Multiple Tenants</a>	• • • • • 24
19.	<a href="#">Case Registration Function ①</a>	• • • • • 25
20.	<a href="#">Case Registration Function ②</a>	• • • • • 26
21.	<a href="#">Case Registration Function ③</a>	• • • • • 27
22.	<a href="#">Case Registration Function ④</a>	• • • • • 28
23.	<a href="#">Stroke Protocol</a>	• • • • • 29
24.	<a href="#">Voice/Video Call Function ①</a>	• • • • • 33
25.	<a href="#">Voice/Video Call Function ②</a>	• • • • • 34
26.	<a href="#">Voice/Video Call Function ③</a>	• • • • • 35
27.	<a href="#">Voice/Video Call Function ④</a>	• • • • • 36
28.	<a href="#">Voice/Video Call Function ⑤</a>	• • • • • 37
29.	<a href="#">Voice Message Function ①</a>	• • • • • 38
30.	<a href="#">Voice Message Function ②</a>	• • • • • 39
31.	<a href="#">Voice Message Function ③</a>	• • • • • 40
32.	<a href="#">Sending Videos ①</a>	• • • • • 41
33.	<a href="#">Sending Videos ②</a>	• • • • • 42
34.	<a href="#">Sending Videos ③</a>	• • • • • 43
35.	<a href="#">Sending Videos ④</a>	• • • • • 44
36.	<a href="#">Sharing PDF</a>	• • • • • 45
37.	<a href="#">Sharing Images</a>	• • • • • 47

# Installation & Registration Process

## ① Install Join from Google Play



Google Play



\*For those who cannot read this QR code please search for “Join allm” on Google Play and install the app that has the icon you can see on the left.

\*Only compatible with Android OS 4.5, or above. If the OS on your device is lower than these please update it.

## ② Start the app

## ③ Tap “Create Account”

Login

Email

Password

Region  
North America

LOGIN

[Forgot password?](#)

CREATE ACCOUNT

## Create Account

### Personal Information

Name

Email

Password



\* 8-20, one number and one letter required.

Repeat Password



### Tenant Information

Tenant Code

Tenant Password



Region

North America



☐ I accept the [Terms of Use](#)

CREATE ACCOUNT

④ Enter your name, email address, and password

\*Symbols can not be used  
Except: "&", "-", "\_", "."

⑤ Enter the tenant code and tenant password

\*If your tenant code or password is unclear, please contact a manager of your tenant.

⑥ Agree to the Terms of Use and tap the checkbox

⑦ Tap "Create Account"

⑧ Enter the code that you received via Email and tap “Verify”

\*If you don't receive the verification code, please tap "Resend the Code" to send again. If you haven't received even after re-sending, please return to the previous screen and confirm the email address you entered is correct.

⑨ Set a 4-digit passcode to unlock the app. And tap “Save”

Setting completed.  
You will be logged in to Join.

## Email Verification

A verification code has been sent to Join\_manual@allm.net

Please enter the verification code to finish creating your account.

VERIFY

RESEND THE CODE

## Passcode Setting

\* 4 digit number.

Your passcode is required to unlock the application after it has not been used for a certain length of time. Please store it safely.

Should you forget your passcode, you will have to reinstall the application, and all of your data will be lost.

Allm Inc. and its affiliates shall not be held responsible or liable in any way should this occur.


SAVE

# Cautions for Using

## ▼ Setting Passcode

For security purposes, if Join hasn't been used for a specified amount of time, it will be locked. A passcode is necessary to unlock it.

① Tap the setting icon



## ■ Change Passcode

### Change PassCode

Please set the passcode by 4 number digits.

Current PassCode

New PassCode

New PassCode(Retry)

CANCEL UPDATE

## ■ Lock Timer (Change time until passcode lock)

### Lock Timer

15 minutes

CANCEL UPDATE

\*If you forget the passcode **you will have to reinstall Join.**

\*Please log in using email address and password that you registered at [Email Login Setup] when you reinstalled the app.

▼ Contact: If you have any questions or concerns please feel free to get in contact with us.

Allm Inc. Join Support Center [support@join-me.biz](mailto:support@join-me.biz)

# Change registered information

## ▼ Change Email

You can change Email address you've registered.




① Tap "Change Email"

② Enter current password and new email address

③ Tap "Change"

④ Enter the code you received

⑤ Tap "Verify"



Email address will be changed.

\*The verification code is valid for 10 mins.



## ▼ Change Password

You can change login password.

① Tap "Change Password"

② Enter current password and new password

③ Tap "CHANGE"



Password will be changed.

※ Password needs to be more than 8 and less than 20 alphanumeric characters.



# Join Reinstallation

## ▼ Email Login

Please do not register with phone number authentication when you reinstalled Join.

Please log in by using **email address and password that you registered at “Account Setting”**.

The screenshot shows the Join app's login interface. At the top, there's a red header with a white pulse icon. Below it, a dashed yellow box highlights the 'Email' and 'Password' input fields. A yellow box highlights the 'Region' dropdown menu, which is currently set to 'North America'. Below these fields is a red 'LOGIN' button. Underneath the button is a link that says 'Forgot password?'. At the bottom of the screen is a white button with a red border that says 'CREATE ACCOUNT'.

① Enter email address & password that you have registered

② Make sure your region is correct

③ Tap “Login”

\*If region setting is incorrect, please set it to the right one by yourself.

\*If you register again, **new account will be created and you can't restore your group settings and talk logs.**

\*When you login using email, your contact settings and group talks will be stored.

\*One-to-one direct chats cannot be restored.

# When you forgot password

## ▼ Reset password

If you forgot password, you can reset new one.

① Tap "Forgot password?"

② Enter Email address and tap "SEND"

③ Access to the URL in the email you receive

④ Enter new password and tap "Apply" to set

**Login**

Email

Password

Region  
North America

LOGIN

Forgot password?

CREATE ACCOUNT

**Reset Password**

To reset your password, please submit your email address below.

Email

Region  
North America

SEND

**Reset Password**

An email has been sent to your address. Please follow the instructions in the email and proceed with resetting your password.

OK

**Join**

**Password reset**

Please set the password again.

**Password**

Please input between 8 and 20 characters

**Re-enter password**


Please input the password again

Apply

\*After setting new password, please go back to the app and proceed with login.

# Contacts and Groups

## ▼ Add Contacts



# Contacts and Groups


## ▼ Create Groups

① Tap "CONTACTS"

② Tap "Create Group"

③ Input group name and tap "Create"


④ Tap the group you've created

⑤ Tap 

⑥ Choose a contact and tap "INVITE"

You can set a group avatar.

Groups are consisted of invited contacts.  
You can communicate with multiple members simultaneously in "group chat"




The screenshots illustrate the process of creating a group:

- Step 1:** The user is on the home screen. The 'CONTACTS' tab is selected.
- Step 2:** A menu is displayed with options: 'Add Contact', 'Create Group', and 'Add Chat Room'. The 'Create Group' option is selected.
- Step 3:** The 'Create Group' screen shows a text input field for the group name and a 'Create' button at the bottom.
- Step 4:** The user is back on the home screen. Under the 'Groups (1)' section, 'Group A' is listed.
- Step 5:** The user taps the 'Add contact' icon (a person with a plus sign) to invite members to the group.
- Step 6:** A 'Select Contacts' modal is shown, allowing the user to choose contacts to invite. 'Dr. smith' is selected, and the 'INVITE' button is tapped.


# Direct Chat

## ▼ Start a direct chat

① Tap  
“CONTACTS”




② Tap a contact you want  
to send a message



③ Tap



④ Input a message and  
tap “Send” to post



Left: Voice Talk  
Right: Video Talk



: Start one-on-one direct chat



: Start voice call



: Start video call

\*You can start voice/video call by tapping the call button in direct chat rooms.

# Direct Chat Functions ①

## ▼ Top Menu



Voice call  
button

Video call  
button


Leave Chat

Notification



You can leave the chat group

You can turn off Notification





When you leave the  
chat will be disappeared  
from the chat list.

# Direct Chat Functions ②

## ▼ Bottom Menu

You can use various functions from the bottom menu.



If you tap stickers in the lowest column, you can see different kind of stickers.

The screenshot shows a form titled "Select DICOM Image" with a red header. The form contains the following fields and options:

- Tenant: Free trial
- Exam. Date: Select Examine Date
- Pt. ID: Input Patient ID
- Pt. Name: Input Patient name
- Type: None
- Rept. No: Input Reception Number

At the bottom of the form are two red buttons: "Reset" and "Search".

You can select a DICOM image and post it to the current chat room. (Explained in the next page)

# Direct Chat Functions ③

## ▼ View & send DICOM images

① Tap "Search"

← Select DICOM

Tenant:  
Free trial

Exam. Date:  
Select Examine Date

Pt. ID:  
Input Patient ID

Pt. Name:  
Input Patient name

Type:  
None

Rcpt. No:  
Input Reception Number

Reset Search

② Choose an examination and tap it

← Select DICOM Image

Exam. Date: Jan 21, 2016	Rcpt. No: 1413804-200221	
PL ID: 1413804-2	Type: MR	Region: HEAD
Pt. Name:		

Exam. Date: Dec 13, 2015	Rcpt. No: 1578501-700221	
PL ID: 1578501-7	Type: CT	Region: HEAD
Pt. Name:		

Exam. Date: Nov 13, 2015	Rcpt. No: 15111301	
PL ID: 15111301	Type: XA	Region: HEAD
Pt. Name:		

Exam. Date: Nov 4, 2015	Rcpt. No: 1930753-500061	
PL ID: 1930753-5	Type: CT	Region: ABDOMEN
Pt. Name:		

Exam. Date: Feb 21, 2011	Rcpt. No: t00001	
PL ID: **001	Type: CR	Region: CHEST
Pt. Name:		

③ Tap "Send"

← Select DICOM Image

Exam. Date:  
Jan 21, 2016

PL ID:  
1413804-2

Pt. Name:

Type:  
MR

Region:  
HEAD

Sex:  
F

Age:

No. of Images:  
497

Time of Exams:  
17:14

Send Preview

④ Tap the thumbnail to view

← Hospital A ER

ok. 20:29

Hi, I am sending a DICOM image to a patient.

Can you have a look? 20:30

Read by 1 20:30

OK, I will check this and get back to you ASAP.

Wed, Jun 22

10:07

Write a message

Preview Screen



The image will be sent on the chat room

\*Users who belong to the same tenant can share and view images in a chat room.

\*You won't be able to see images of different tenants unless managers set permission.


\*An error message will be displayed if you don't have relevant permission.




# Group Chat

## ▼ Start a group chat


① Tap a group name





② Tap



③ Input a text message and tap "Send" to post



If you tap the aviator...




If the user is on your contact list, user detail screen will be shown.

If not, you can add him/her by tapping "Add Contact".

# Group Chat Functions ①

## ▼ Top Menu



You can change group name and avatar, invite and delete members, etc.

Leave the group

It will be displayed only if a group site is set by management function.


Tap to play



Web page will be displayed on a browser if a link to the page is set.


You can invite your contacts to a group.

Tick a contact and tap "Select"



You can turn on/off notification of each group.

Groups you turned off notification can be checked on chats list.



# Group Chat Functions ②

## ▼ Bottom Menu

You can use various functions from the bottom menu.

The diagram illustrates the bottom menu of a group chat interface, showing various functions available to users. The central element is a chat input area with a "Write a message" placeholder and a microphone icon. Surrounding this are several function callouts:

- Input text messages**: Points to the text input field.
- Select photos/videos**: Points to the gallery icon.
- Take photos/videos**: Points to the camera icon.
- Stroke Protocol**: Points to the protocol icon.
- Display stickers**: Points to the sticker icon.
- Medical Information**: Points to the medical information icon.
- Record voice messages**: Points to the microphone icon.

Below the main chat area, three additional screens are shown:

- Sticker Selection Screen**: Displays a grid of cartoon cow stickers with various expressions and text like "ON DUTY!", "SO TIRED", "DON'T FORGET TO TAKE YOUR BREAK", "OK!!", "FIGHT", and "I'M SLEEPING".
- Select DICOM Image Screen**: A form for selecting a DICOM image. It includes fields for Tenant (Free trial), Exam Date (Select Examine Date), Pt. ID (Input Patient ID), Pt. Name (Input Patient name), Type (None), Rpt. No., and Input Reception Number. There are "Reset" and "Search" buttons at the bottom.
- Stroke Protocol Screen**: A form for recording a case event. It includes fields for ER Type (Stroke), Patient, Protocol, and DateTime (Jun 22, 2017 6:23 PM). There is a "Protocol" button at the bottom.

Two callouts provide additional context:

- You can quickly post an case event. (Explained later)**: Points to the "Protocol" button on the Stroke Protocol screen.
- If you tap "DICOM", you can select a DICOM image and post it to the current chat room. (Explained in the next page)**: Points to the "DICOM" button on the Select DICOM Image screen.

# Group Chat Functions ③

## ▼ View & send DICOM images

① Tap "Search"

② Choose an examination and tap it

③ Tap "Send"

④ Tap the thumbnail to view

Preview Screen

← Select DICOM Image

Tenant:  
Free trial

Exam. Date:  
Select Examination Date

PL ID:  
Input Patient ID

PL Name:  
Input Patient name

Type:  
None

Rcpt. No:  
Input Reception Number

Reset Search

← Select DICOM Image

Exam. Date: Jan 21, 2016	Rcpt. No: 1413804-200221	
PL ID: 1413804-2	Type: MR	Region: HEAD
PL Name:		

Exam. Date: Dec 13, 2015	Rcpt. No: 1578501-700221	
PL ID: 1578501-7	Type: CT	Region: HEAD
PL Name:		

Exam. Date: Nov 13, 2015	Rcpt. No:	
PL ID: 15111301	Type: XA	Region: ABDOMEN
PL Name:		

Exam. Date: Nov 4, 2015	Rcpt. No: 1930753-500061	
PL ID: 1930753-5	Type: CT	Region: ABDOMEN
PL Name:		

Exam. Date: Feb 21, 2011	Rcpt. No: t00001	
PL ID: **001	Type: CR	Region: CHEST
PL Name:		

← Select DICOM Image

Exam. Date: Jan 21, 2016	Rcpt. No: 1413804-200221
PL ID: 1413804-2	DOB:
PL Name:	

Type: MR	Sex: F	Age:
No. of Images: 497	Time of Exams: 17:14	

Send Preview

← Group A

18:13 アルム病院 アルム太郎 j... group.

18:15 Dr. A removed User A fr... group.

18:16 User A joined to th...

User A

Test message 18:16

Thu, Aug 24

10:55

Write a message



\*Users who belong to the same tenant can share and view images in a chat room.

\*You won't be able to see images of different tenants unless managers set permission.

\*An error message will be displayed if you don't have relevant permission.

The image will be sent on the chat room

# Notifications①

## ▼ Turn off notifications


You can turn off notifications not to receive any notifications from Join.

① Tap the setting icon

② Tap "Notifications"

③ Tap "Allow Notifications"

OFF




# Notifications②


## ▼ Change notification sound

You can select and set a notification sound you like from the list.


① Tap "Notifications"



② Tap "Notification Sound"



③ Select a sound and tap the icon




The selected sound will be set

\*The selected notification sound will be applied to any notification sent from Join.

\*The sound can not be set per group.

# DICOM Image Acquisition

## ▼ DICOM viewer



\*Users who belong to the same tenant can share and view images in a chat room.

\*You won't be able to see images of different tenants unless managers set permission.

\*An error message will be displayed if you don't have relevant permission.

# Belonging to Multiple Tenants

## ▼ Register multiple tenants

① Tap Settings or your display name

② Tap “Add Tenant”

③ Enter tenant code and tap “Send Application”

Tenant request will be sent to the administrator

After the administrator approved, you can join the tenant


\*It may take a while to be reflected after sending a request.  
You need to reboot the app to update display of your tenants in Settings.




# Case Registration Functions ①

## ▼ Register an case event

① Tap “CASE LIST”




② Tap “+” to register




③ Enter information

④ Select a group and tap “OK”



⑤ Tap “DONE” to save




\*If you don't have permission to add cases "+" button will not be displayed.

Registered information will be posted to the group chat when you save it.



# Case Registration Functions ②

## ▼ View patient information



# Case Registration Functions ③

## ▼ Start and stop tracking



\*Tracking is set to off by default.  
You need to enable tracking on the Settings.

# Case Registration Functions ④

## ▼ Edit patient information

① Tap "Edit"

② Edit patient information

③ Tap "DONE"

The first screenshot shows the 'Patient Info' screen with a red header and a red edit icon in the top right corner. The patient information includes: Patient ID: Test0001, Sex: Male, Age: 85, Weight: 58.20, Patient Type: Transfer, NIHSS: 21, ER Type: Stroke, Group: グループ A, Onset/LKW: , and Note: . At the bottom are buttons for 'Start Tracking', 'View Location', 'Time Stamps', and 'DICOM'.

The second screenshot shows the 'Edit' screen with a red header and a red 'DONE' button in the top right corner. The patient information is the same as the first screenshot, but the 'Onset / LKW' field is now populated with 'Aug 24, 2017 14:50' and the 'Note' field contains 'xxxx \*\*\* Blvd, \*\*\* City'.


The third screenshot shows a group chat interface with a red header and a red edit icon in the top right corner. The chat history shows a message from 'Registered by Dr. A' with the patient information. A new message from 'Updated by Dr. A' is being composed, containing the updated patient information: Tenant: Free trial, Patient ID: allmtest0824, ER Type: Stroke, Sex: Male, Weight: 50, Age: 65, Patient Type: Transfer, On Set / LKW: Aug 24, 2017 2:50 PM, NIHSS: 20, and Note: xxxx \*\*\* Blvd, \*\*\* City. The chat interface includes a 'Write a message' field and various icons for attachments and actions.

Edited items are automatically posted to the group chat when you save it.


# Protocol

## ▼ Post an case event

① Tap




② Tap “Patient” and choose a patient



Tenant	Last Update
スキルアップジャパン株式会社	Feb 19, 2018 17:56
Patient ID: test0002	Sex: Male, Age: 65
Tenant: Free trial	Last Update: Oct 5, 2017 17:53
Patient ID: allmtest0824	Sex: Male, Age: 65
Tenant: スキルアップジャパン株式会社	Last Update: Aug 9, 2017 16:36
Patient ID: test0809	Sex: Male, Age: 50
Tenant: スキルアップジャパン株式会社	Last Update: Jun 22, 2017 18:53
Patient ID: AllmTest001	Sex: Male, Age: 65
Tenant: スキルアップジャパン株式会社	Last Update: Jun 22, 2017 18:47
Patient ID: AllmTest002	Sex: Female, Age: 52
Tenant: スキルアップジャパン株式会社	Last Update: Jun 21, 2017 12:35

③ Tap “Protocol”




ER Type: Stroke

Protocol: Select protocol

DateTime: Feb 19, 2018 18:36

④ Tap a corresponded protocol



Select protocol

- tPA
- Arrival to hospital
- CT/MRI start
- CT/MRI end
- Into Angio Suite
- Puncture
- Lesion / First pass

\*Protocol is available only in group chats.

\*Protocol button won't be displayed in direct chats.

\*You can select a corresponding protocol to the case

\*Only cases registered to the group or transferred from another group can be selected.

⑤ Tap date time

Protocol

Done

Patient

Tenancy	Fre	Last Update	Oct 5, 2017 17:53
Patient	allm	Sex	Male
		Age	65

ER Type

Stroke

Protocol

Arrival to hospital

DateTime

Jun 22, 2017 18:23

⑥ Set event occurrence time and tap "OK"

Emergency case event

2017

Thu, Jun 22

June 2017

S	M	T	W	T	F	S
				1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	

CANCEL OK

6:41 AM

CANCEL

Drag the red circle to select time

Tap to revise

6:41 AM

CANCEL OK

⑦ Tap "DONE"

Protocol

Patient

Tenant	Last Update
Free trial	Oct 5, 2017 17:53

Patient ID	Sex	Age
allmtest0824	Male	65

ER Type

Stroke

Protocol

Arrival to hospital

DateTime

Jun 22, 2017 18:41

DONE

Group A

Age: 65  
Patient Type: Transfer

Registered by User B

18:36


Tenant: スキルアップジャパン株式会社  
Patient ID: AllmTest002  
ER Type: Stroke  
Sex: Female  
Weight: 48  
Age: 52  
Patient Type: Transfer  
On Set / LKW: Jun 22, 2017 18:36

Registered by User B


18:44  
Arrival to hospital  
AllmTest002  
Jun 22, 2017 18:41  
Registered by User B

Write a message

Notification  
on iOS



Notification  
on Android




Event, Patient ID,  
and date time are  
posted to the chat


When a Protocol is posted,  
all members of the same group will receive  
a pop-up notification

## ▼ View protocols


① Tap patient information




Case List



② Tap "Time Stamps"



You can add Time Stamps by tapping 「+」 here




Time Stamps posted on the group chats are listed for each patient.



# Voice/Video Call Functions ①

## ▼ Start voice/video call

【Contact list】




# Voice/Video Call Functions ②

## ▼ Start voice/video call

### 【Direct chat】


Tap the call button to start voice/video call



### 【Calls tab】


① Tap calls tab to display call history

② Tap the phone button to start call



# Voice/Video Call Functions ③

## ▼ In coming call screen



You can reject the call.  
A message which says "User is busy at the moment." will be displayed on caller's screen.

# Voice/Video Call Functions ④


## ▼ Voice call screen

Voice call screen

Switch display

You can switch display during a call.

Caller's Display name and the call duration




Function buttons

Turn on/off microphone

Turn on/off camera

Turn on/off speaker



Hang up button

# Voice/Video Call Functions ⑤

## ▼ Video call screen

Video call screen



Function buttons



Switch front/back cameras



Switch display



Turn on/off microphone

Hang up button

Turn on/off camera

Your face on screen

# Voice Message Functions ①

## ▼ Record and post a voice message

Direct / Group Chat

① Tap

② Hold

③ Speak while holding the button.

④ Leave your finger from the button.

Good morning!

A voice message is posted on the chat room.


The image illustrates the process of recording and posting a voice message in a chat room through four sequential steps:

- ① Tap:** The user taps the microphone icon at the bottom right of the chat interface.
- ② Hold:** The user holds the red circular microphone button.
- ③ Speak while holding the button:** The user speaks while holding the button, and the recorded voice is shown as a blue scribble.
- ④ Leave your finger from the button:** The user releases the button, and the recorded voice message "Good morning!" is posted to the chat.


The chat interface shows a group chat for "Group A" with various medical notes and a "Read by 1" status.

# Voice Message Functions ②

## ▼ Cancel recording




① Keep holding the button and move your finger outside the button.



The outer line turns red.


② Make sure the outer line turned red and leave your finger.



It stops recording and a voice message won't be sent.

# Voice Message Functions ③

## ▼ Play a voice message






# Sending Videos ①

## ▼ Record and post a video

Direct / Group Chat

① Tap 

② Tap "Record a video"

③ Tap the Record button to start recording.

④ Check the video and tap "✓"

\*Display will be different depends on a device you use.


Confirmation screen will be displayed if you stop recording or 5 minutes have passed.

The video is sent and a thumbnail with a play icon (▶) is displayed.

# Sending Videos ②

## ▼ Upload a video file

Direct / Group Chat



\*Videos longer than 5 minutes can not be uploaded.

Please trim the video with another app to make it **shorter than 5 minutes** and send it.

\*Please read User Manual for your device to see how to edit videos.

\*You might need to download an editing app depending on the device.


The video is sent and a thumbnail with a play icon(▶) is displayed.

# Sending Videos ③

## ▼ Play a video

Direct / Group Chat


① Tap a thumbnail with



② Tap the screen to change the display of user interface.

On

Off



The video is played automatically.

③ Tap the back key to go back to the chat room.




Play/Pause


Control Volume

# Sending Videos ④

## ▼ Adjust quality of Photos and Videos



### Quality Comparison



\*If you enable HD upload, the quality rises and the data file grows. It takes longer to send/receive videos/photos.

# Sharing PDF

## ▼ Share PDF files


You can share PDF file on the Internet or the file you've received through other applications to Join.

① Open the file

Web browser




Other application



② Tap "Send file"




③ Tap Join




\*Availability of this feature depends on applications you share files from.

④ Select a chatroom and tap



⑤ Check the file and tap "SEND"



The PDF file will be shared in the chatroom.

# Sharing images

## ▼ Share images

You can share images from other applications to Join.

① Open the file

② Tap "Share" or  
 in menu




③ Tap Join

E.g.

Gmail

Gallery



Web browser



\*Any image format will be converted to JPEG and sent to Join.  
GIF images will be shared as non-animated images.

\*Availability of this feature depends on applications you share files from.

④ Select a chatroom and tap



The image will be shared in the chatroom.